Human GIP Protein


Cat. No. GIP-HM201

ϗͶϲϿ·Ͷϩ

Description	
Source	Recombinant Human GIP Protein is expressed from HEK293 with hFc tag at the C-Terminus.
	It contains Glu22-Gln93.
Accession	P09681
Molecular Weight	The protein has a predicted MW of 34.9 kDa. Due to glycosylation, the protein migrates to 40-50 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC
Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Background	
	The potential application of glucose-dependent insulinotropic polypeptide (gastric inhibitory polypeptide, GIP) in the management of obesity and type 2 diabetes has been controversial. Initial interest in the therapeutic use of GIP was dampened by evidence that its insulinotropic activity was reduced in type 2 diabetes and by reports that it increased glucagon secretion and adipose deposition in non-diabetic individuals.


Assay Data

Tris-Bis PAGE



Human GIP on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



1



Human GIP Protein

KVCJUS

The purity of Human GIP is greater than 95% as



Immobilized Human GIP, hFc Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Anti-GIP Antibody, hFc Tag with the EC50 of 88.1ng/ml determined by ELISA.

0.0-