

Mouse Osteoactivin/GPNMB Protein

Cat. No. GNB-MM101

Description

Source	Recombinant Mouse Osteoactivin/GPNMB Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Lys23-Asn502.
Accession	NP_444340.3
Molecular Weight	The protein has a predicted MW of 54.7 kDa. Due to glycosylation, the protein migrates to 80-110 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per ug by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

GPNMB is highly expressed in macrophages and microglia, which are cells involved with innate immune response in the periphery and the brain, respectively. Some studies have shown increased levels of GPNMB in pro-inflammatory conditions, such as LPS treatment, and in pathological conditions, such as neurodegenerative diseases and cancer.

Assay Data


Bis-Tris PAGE



Mouse Osteoactivin on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse Osteoactivin is greater than 95% as determined by SEC-HPLC.