

FITC-Labeled Human Her3/ErbB3 Protein

Cat. No. HER-HM403F

Description

Source	Recombinant FITC-Labeled Human Her3/ErbB3 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Ser20-Thr643.
Accession	P21860-1
Molecular Weight	The protein has a predicted MW of 71.6 kDa. Due to glycosylation, the protein migrates to 72-75 kDa based on Bis-Tris PAGE result.
Wavelength	Excitation Wavelength: 490 nm Emission Wavelength: 520 nm
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 6 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Her3, also called ErbB3, is a type I membrane glycoprotein that is a member of the ErbB family of tyrosine kinase receptors. Her3 is expressed in keratinocytes, melanocytes, skeletal muscle cells, embryonic myoblasts and Schwann cells. Monomeric Her3 serves as a low affinity receptor for the heregulins (HRG).

Assay Data

Bis-Tris PAGE




FITC-Labeled Human Her3 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

FITC-Labeled Human Her3, His Tag ELISA

0.5µg FITC-Labeled Human Her3, His Tag Per Well



Immobilized FITC-Labeled Human Her3, His Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Human NRG1 Beta 1, hFc Tag with the EC50 of 23.6ng/ml determined by ELISA.