

Biotinylated Human/Cynomolgus Her4/ErbB4 Protein, Ultra Low Endotoxin

Cat. No. HER-HM4B4B-UL

Description	
Source	Recombinant Biotinylated Human/Cynomolgus Her4/ErbB4 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Gln26-Pro651.
Accession	Q15303-1(Human) / A0A2K5W0D7(Cynomolgus)
Molecular Weight	The protein has a predicted MW of 72.7 kDa. Due to glycosylation, the protein migrates to 75-110 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Her4, a member of the EGF receptor family, plays a variety of roles in physiological and pathological states. Genetic studies have indicated a link between Her4 and type 2 diabetes and obesity. Her4 may play an important role in glucose homeostasis and lipogenesis. Her4 deficiency-related obesity and adipose tissue inflammation may contribute to the development of MetS.

Assay Data

Bis-Tris PAGE




Biotinylated Human/Cynomolgus Her4 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

Assay Data

Biotinylated Human/Cynomolgus Her4, His Avi Tag ELISA

0.2µg Biotinylated Human/Cynomolgus Her4, His Avi Tag Per Well



Immobilized Biotinylated Human/Cynomolgus Her4, His Tag at 2µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Human NRG1 Beta 1, hFc Tag with the EC50 of 8.8ng/ml determined by ELISA.