

Mouse IgG2C Protein, Ultra Low Endotoxin

Cat. No. IGG-MM02C-UL

Description

Source	Recombinant Mouse IgG2C Protein is expressed from HEK293 without tag. It contains Ile97-Ser330.
Accession	F6TQW2
Molecular Weight	The protein has a predicted MW of 26.10 kDa. Due to glycosylation, the protein migrates to 30-40 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Mouse IgG2c is a subclass of mouse antibody IgG, it can bind FcγRIV with high affinity and the Fc loop residues of mouse IgG2c promote greater receptor-binding affinity than mouse IgG2b or human IgG1.

Assay Data


Bis-Tris PAGE



Mouse IgG2C on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse IgG2C is greater than 95% as determined by SEC-HPLC.