

Mouse IL-1R2/IL-1 RII/CD121b Protein

Cat. No. IL1-MM1R2

Description

Source	Recombinant Human IL-1R2/IL-1 RII/CD121b Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Phe14-Glu355.
Accession	P27931
Molecular Weight	The protein has a predicted MW of 39.2 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage

Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3-6 months after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Interleukin (IL)-1 has been reported to be involved in the development of tuberculosis (TB). IL1R1 and IL1R2 encode a cytokine receptor that belongs to the IL-1 receptor family. IL1R2 (Interleukin 1 Receptor Type 2) is a Protein Coding gene. Diseases associated with IL1R2 include Endometriosis and Mastitis. Among its related pathways are Hematopoietic cell lineage and IL-1 signaling pathway. Gene Ontology (GO) annotations related to this gene include interleukin-1 receptor activity and interleukin-1, type II, blocking receptor activity.

Assay Data


Bis-Tris PAGE



Mouse IL-1R2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Mouse IL-1R2 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Mouse IL-1R2, His Tag ELISA

0.5µg Mouse IL-1R2, His Tag Per Well



Immobilized Mouse IL-1R2, His Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Mouse IL-1 Beta, His Tag with the EC50 of 1.2µg/ml determined by ELISA (QC Test).