

Biotinylated Cynomolgus IL-25/IL-17E Protein

Cat. No. IL2-CM425B

Description

Source	Recombinant Biotinylated Cynomolgus IL-25/IL-17E Protein is expressed from HEK293 with His tag at the C-terminus and Avi tag at the N-terminus. It contains Tyr33-Gly176.
Accession	G7P9V9
Molecular Weight	The protein has a predicted MW of 19.99 kDa. Due to glycosylation, the protein migrates to 25-32 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 µm filtered solution in PBS, 8% trehalose (pH 7.4).
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

IL-25, also known as IL-17E, is a member of the IL-17 cytokine family mostly produced by epithelial cells and innate immune cells. After binding to the IL-17RB/IL-17RA complex, IL-25 induces downstream signaling responses in epithelial cells and type 2 lymphocytes, which initiates, propagates, and sustains type 2 immunity. The function of IL-25 in allergic diseases such as asthma has been well established, and now also is extended to diseases such as inflammatory bowel disease and cancer.


Assay Data

Bis-Tris PAGE



Biotinylated Cynomolgus IL-25 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC



The purity of Biotinylated Cynomolgus IL-25 is greater than 95% as determined by SEC-HPLC.

Assay Data

SEC-MALS




The purity of Biotinylated Cynomolgus IL-25 is greater than 95% and the molecular weight of this protein is around 40-65 kDa as determined by SEC-MALS.

ELISA Data

Biotinylated Cynomolgus IL-25, His Tag ELISA

0.5µg Human IL-17RB, His Tag Per Well




Immobilized Human IL-17RB, His Tag (Cat. IL1-HM1RB) at 5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Cynomolgus IL-25, His Tag with the EC50 of 11.1ng/ml determined by ELISA.

Log Biotinylated Cynomolgus IL-25, His Tag Conc.(µg/ml)

ELISA Data

Biotinylated Cynomolgus IL-25, His Tag ELISA

0.05µg Biotinylated Cynomolgus IL-25, His Tag Per Well



Immobilized Biotinylated Cynomolgus IL-25, His Tag at 0.5µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Anti IL-25 Antibody, hFc Tag with the EC50 of 4.4ng/ml determined by ELISA.

Log Anti IL-25 Antibody, hFc Tag Conc.(µg/ml)