

# Biotinylated Human IL-23 alpha&IL-12 beta Protein

Cat. No. IL2-HM4ABB

## Description

<b>Source</b>	Recombinant Biotinylated Human IL-23 alpha&IL-12 beta Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Arg20-Pro189 (IL-23 alpha) and Ile23-Ser328 (IL-12 beta).
<b>Accession</b>	Q9NPF7(IL-23 alpha)&P29460(IL-12 beta)
<b>Molecular Weight</b>	The protein has a predicted MW of 21.7 kDa (IL-23 alpha) and 34.7 kDa (IL-12 beta). Due to glycosylation, the protein migrates to 24 kDa (IL-23 alpha) and 40-50 kDa (IL-12 beta) based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 1EU per µg by the LAL method.
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage


<b>Formulation</b>	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
<b>Storage</b>	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

The interleukin (IL)-12 family of cytokines, including IL12 and IL 23, play an important role in driving aberrant Th1 and Th17 immune responses in patients with Crohn's disease (CD). Targeting this pathway has opened new avenues for therapeutic intervention. The balance of proinflammatory cytokines interleukin (IL)-12 and IL-23 plays a key role in shaping the development of antitumor or protumor immunity.

## Assay Data

### Bis-Tris PAGE




Biotinylated Human IL-23 alpha&IL-12 beta on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC

# Biotinylated Human IL-23 alpha&IL-12 beta Protein

Cat. No. IL2-HM4ABB

## Assay Data




The purity of Biotinylated Human IL-23 alpha&IL-12 beta is greater than 95% as determined by SEC-HPLC.

## ELISA Data

### Biotinylated Human IL-23 alpha&IL-12 beta, His Tag ELISA

0.2µg Biotinylated Human IL-23 alpha&IL-12 beta, His Tag Per Well




Immobilized Biotinylated Human IL-23 alpha&IL-12 beta, His Tag at 2µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Human IL-23R, hFc Tag with the EC50 of 21.0ng/ml determined by ELISA (QC Test).

## ELISA Data

### Biotinylated Human IL-23 alpha&IL-12 beta, His Tag ELISA

0.05µg Biotinylated Human IL-23 alpha&IL-12 beta, His Tag Per Well



Immobilized Biotinylated Human IL-23 alpha&IL-12 beta, His Tag at 0.5µg/ml (100µl/Well) on streptavidin (5µg/ml) precoated plate. Dose response curve for Anti-IL-23 Antibody, hFc Tag with the EC50 of 2.3ng/ml determined by ELISA.