

Biotinylated Mouse IL-2 Protein

Cat. No. IL2-MM401B

Description

Source	Recombinant Biotinylated Mouse IL-2 Protein is expressed from HEK293 with His tag and Avi tag at the N-Terminus. It contains Ala21-Gln169.
Accession	P04351-1
Molecular Weight	The protein has a predicted MW of 20.1 kDa. Due to glycosylation, the protein migrates to 27-30 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per ug by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Interleukin-2 (IL-2) is an interleukin, a type of cytokine signaling molecule in the immune system. It is a 15,5 - 16 kDa protein that regulates the activities of white blood cells (leukocytes, often lymphocytes) that are responsible for immunity. IL-2 has essential roles in key functions of the immune system, tolerance and immunity, primarily via its direct effects on T cells.

Assay Data


Bis-Tris PAGE



Biotinylated Mouse IL-2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

Biotinylated Mouse IL-2, His-Avi Tag ELISA
0.5µg Mouse IL-2 R alpha, His Tag Per Well



Immobilized Mouse IL-2 R alpha, His Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Mouse IL-2, His Tag with the EC50 of 0.15µg/ml determined by ELISA (QC Test).