Mouse IL-15RA/IL-15 R alpha/CD215 Protein

Cat. No. IL5-MM2RA



Description	
Source	Recombinant Mouse IL-15RA/IL-15 R alpha/CD215 Protein is expressed from HEK293 with hFc tag at the C-Terminus.
	It contains Gly33-Lys205.
Accession	Q60819-1
Molecular Weight	The protein has a predicted MW of 45.1 kDa. Due to glycosylation, the protein migrates to 75-82 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC
Formulation and	Storage
	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trebalose is added as protectant before


Formulation	lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Interleukin-15 receptor alpha (IL-15R alpha) is a high affinity IL-15 binding protein that is crucial for mediating IL-15 functions such as memory CD8 T cell proliferation and NK, NK/T cell, and intestinal intraepithelial lymphocyte development.

Assay Data

Tris-Bis PAGE




Mouse IL-15RA on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

KAGTUS


Assay Data



The purity of Mouse IL-15RA is greater than 95% as determined by SEC-HPLC.

ELISA Data

Mouse IL-15RA, hFc Tag ELISA 0.1µg Human IL-15, No Tag Per Well



Immobilized Human IL-15, No Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Mouse IL-15RA, hFc Tag with the EC50 of 63.4ng/ml determined by ELISA.