

Human IL-5 R alpha/CD125 Protein, Ultra Low Endotoxin



Cat. No. ILR-HM45R-UL

Description	
Source	Recombinant Human IL-5 R alpha/CD125 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Asp21-Glu335.
Accession	Q01344-1
Molecular Weight	The protein has a predicted MW of 38.7 kDa. Due to glycosylation, the protein migrates to 50-68 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Interleukin5 Receptor alpha (IL5 R alpha), also known as CD125, is a 60 kDa hematopoietin receptor that plays a dominant role in eosinophil biology. This is the receptor for interleukin-5. The alpha chain binds to IL5.	

Assay Data

Bis-Tris PAGE



Human IL-5R alpha on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Human IL-5 R alpha/CD125 Protein, Ultra Low Endotoxin

Cat. No. ILR-HM45R-UL




Assay Data



The purity of Human IL-5R alpha is greater than 95% as determined by SEC-HPLC.

SPR Data



Human IL-5, His-Avi Tag immobilized on CM5 Chip can bind Human IL-5 R alpha, His-Avi Tag with an affinity constant of 5.70 nM as determined in SPR assay (Biacore T200).