

Mouse LRG1 Protein, Ultra Low Endotoxin



Cat. No. LRG-MM101-UL

Description	
Source	Recombinant Mouse LRG1 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Leu33-Leu342.
Accession	Q91XL1
Molecular Weight	The protein has a predicted MW of 34.9 kDa. Due to glycosylation, the protein migrates to 48-52 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Diabetic nephropathy (DN) is an important public health concern of increasing proportions and the leading cause of end-stage renal disease (ESRD) in diabetic patients. It is one of the most common long-term microvascular complications of diabetes mellitus that is characterized by proteinuria and glomerular structural changes. LRG1 is a novel pro-angiogenic factors involved in the abnormal angiogenesis and renal fibrosis in DN.

Assay Data


Bis-Tris PAGE



Mouse LRG1 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse LRG1 is greater than 95% as determined by SEC-HPLC.