

# Human LRP10 Protein, Ultra Low Endotoxin

Cat. No. LRP-HM110-UL

Description	
Source	Recombinant Human LRP10 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains His17-Lys440.
Accession	Q7Z4F1-1
Molecular Weight	The protein has a predicted MW of 47.14 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

**Background**

LDL receptor-related protein (LRP) 10 was recently identified as a Parkinson's disease gene through genome-wide linkage and sequencing analysis, but its role in Parkinson's disease in various populations is still unclear.

## Assay Data

### Bis-Tris PAGE



Human LRP10 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.