

Human LRP1 (26-109) Protein

Cat. No. LRP-HM611



Description

Source	Recombinant Human LRP1 (26-109) Protein is expressed from HEK293 with His tag and Strep II tag at the C-terminus.
	It contains Thr26-Arg109.
Accession	Q07954-1
Molecular Weight	The protein has a predicted MW of 11.74 kDa. Due to glycosylation, the protein migrates to 13-20 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage

Formulation	Lyophilized from 0.22 µm filtered solution in 10mM HEPES, 150mM NaCl, 5mM CaCl ₂ , 0.005% Tween-20 (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Low-density lipoprotein (LDL)-related protein-1 (LRP-1) is a member of LDL receptor family that is implicated in lipoprotein metabolism and in the homeostasis of proteases and protease inhibitors. In addition to its function as a scavenger receptor for various ligands, LRP-1 has been shown to transduce multiple intracellular signal pathways including mitogen-activated protein kinase (MAPK), Akt, Rho, and the integrin signaling. LRP-1 signaling plays an important role in the regulation of diverse cellular process, such as cell proliferation, survival, motility, differentiation, and transdifferentiation, and thus participates in the pathogenesis of organ dysfunction and injury.

Assay Data

Bis-Tris PAGE



Human LRP1 (26-109) on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.