

FITC-Labeled Human HLA-A*11:01&B2M&KRAS G12V (VVGAVGVGK) Tetramer Protein



Cat. No. MHC-HE006TF

Description

Source	Recombinant FITC-Labeled Human HLA-A*11:01&B2M&KRAS G12V (VVGAVGVGK) Tetramer Protein is expressed from E.coli with His tag and Avi tag at the C-terminus. FITC-Labeled Human HLA-A*11:01&B2M&KRAS G12V (VVGAVGVGK) Tetramer is assembled by biotinylated monomer and FITC-labeled streptavidin. It contains Gly25-Thr305(HLA-A*11:01), Ile21-Met119(B2M) and VVGAVGVGK peptide.
Accession	AAV53343.1(HLA-A*11:01)&P61769(B2M)&VVGAVGVGK
Molecular Weight	The protein has a predicted MW of 35.36 kDa (HLA-A*11:01) and 11.9 kDa (B2M) same as Bis-Tris PAGE result.
Wavelength	Excitation Wavelength: 490 nm Emission Wavelength: 520 nm
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in 20mM Tris, 200mM NaCl (pH 8.0).
Storage	Valid for 6 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) is the most commonly mutated oncogene in human cancer. The developments of many cancers depend on sustained expression and signaling of KRAS, which makes KRAS a high-priority therapeutic target. The virtual screening approach to discover novel KRAS inhibitors and synthetic lethality interactors of KRAS are discussed in detail.

Assay Data


Bis-Tris PAGE



FITC-Labeled Human HLA-A*11:01&B2M&KRAS G12V (VVGAVGVGK) Tetramer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of FITC-Labeled Human HLA-A*11:01&B2M&KRAS G12V (VVGAVGVGK) Tetramer is greater than 95% as determined by SEC-HPLC.