

Biotinylated Human HLA-A*02:01&B2M&MAGE-A1 (KVLEYVIKV) Monomer Protein



Cat. No. MHC-HE473B

Description

Source	Recombinant Biotinylated Human HLA-A*02:01&B2M&MAGE-A1 (KVLEYVIKV) Monomer Protein is expressed from E.coli with His tag and Avi tag at the C-terminus. It contains Gly25-Thr305(HLA-A*02:01), Ile21-Met119(B2M) and KVLEYVIKV peptide.
Accession	A0A140T913(HLA-A*02:01)&P61769(B2M)&KVLEYVIKV
Molecular Weight	The protein has a predicted MW of 35.6 kDa (HLA-A*02:01) and 11.9 kDa (B2M) same as Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in 20mM Tris, 200mM NaCl (pH 8.0).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

MAGE-A1 belongs to the chromosome X-clustered genes of cancer-testis antigen family and is normally expressed in the human germ line but is also overexpressed in various tumors.


Assay Data

Bis-Tris PAGE



Biotinylated Human HLA-A*02:01&B2M&MAGE-A1 (KVLEYVIKV) Monomer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Biotinylated Human HLA-A*02:01&B2M&MAGE-A1 (KVLEYVIKV) Monomer is greater than 95% as determined by SEC-HPLC.