

Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTRILTI) Monomer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM409B-UL

Description	
Source	Recombinant Biotinylated Human HBV(HLA-A*02:01) Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus It contains Gly25-Thr305(HLA-A*02:01),Ile21-Met119(B2M) and FLLTRILTI peptide.
Accession	A0A140T913(HLA-A*02:01)&P61769(B2M)&FLLTRILTI peptide
Molecular Weight	The protein has a predicted MW of 50.5 kDa. Due to glycosylation, the protein migrates to 52-62 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt.-80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Hepatitis B virus (HBV), is the leading cause of liver diseases infecting an estimated 240 million persons worldwide. The HBV prevalence rates are variables between different countries, with an high level of endemicity in the south-eastern part of Europe. Seven main HBV-D subgenotypes have been described until now (D1-D7).

Assay Data

Bis-Tris PAGE



Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTRILTI) Monomer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTRILT) Monomer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM409B-UL

Assay Data



The purity of Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTRILT) Monomer is greater than 95% as determined by SEC-HPLC.