Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer Protein





Description	
Source	Recombinant Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus.
	It contains Gly25-Thr305(HLA-A*11:01) and Ile21-Met119(B2M).
Accession	AAV53343.1(HLA-A*11:01)&P61769(B2M)
Molecular Weight	The protein has a predicted MW of 48.50 kDa. Due to glycosylation, the protein migrates to 50-60 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC
Formulation and Storage	
Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before


Formulation	lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Peptide Ready HLA-A*11:01&B2M Monomer is absent from peptide, namely peptide-receptive MHC. It can be loaded with antigenic peptides matching HLA-A*11:01. Peptide ready MHC molecules comprising human HLA alleles and B2M, which can be readily tetramerized and loaded with peptides of choice in a high-throughput manner.

Assay Data

Tris-Bis PAGE



SEC-HPLC


Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer Protein

Cat. No. MHC-HM41RB



Assay Data



The purity of Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer is greater than 95% as determined by SEC-HPLC.

Biotinylated Human Peptide Ready HLA-A*11:01&B2M Monomer Protein

Cat. No. MHC-HM41RB




Assay Data

ELISA Data

Biotinlylated Human Peptide Ready HLA-A*11:01&B2M Monomer, His Tag ELISA

0.05μg Biotinlylated Human Peptide Ready HLA-A*11:01&B2M Monomer, His Tag Per Well



Immobilized Biotinlylated Human Peptide Ready HLA-A*11:01&B2M Monomer, His Tag at 0.5µg/ml (100µl/Well) on streptavidin(5µg/ml) precoated plate. Dose response curve for Anti-HLA class I (W6/32) Antibody, hFc Tag with the EC50 of 2.8ng/ml determined by ELISA.