

Human HLA-A*24:02&B2M&Survivin 2B (AYACNTSTL) Monomer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM430-UL

Description

Source	Recombinant Human Survivin 2B(HLA-A*24:02) Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Gly25-Thr305(HLA-A*24:02), Ile21-Met119(B2M) and AYACNTSTL peptide.
Accession	AAA59600.1(HLA-A*24:02)&P61769(B2M)&AYACNTSTL
Molecular Weight	The protein has a predicted MW of 50.20 kDa. Due to glycosylation, the protein migrates to 53-63 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Survivin-2B, a known splice variant of survivin, has been reported to promote cell death in some cancer cells, although it keeps prosurvival function in others. survivin-2B promoted autophagy and further regulated cell death by accumulating and stabilizing IKK alpha in the nucleus.

Assay Data

Bis-Tris PAGE



Human HLA-A*24:02&B2M&Survivin 2B (AYACNTSTL) Monomer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Human HLA-A*24:02&B2M&Survivin 2B (AYACNTSTL) Monomer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM430-UL

Assay Data



The purity of Human HLA-A*24:02&B2M&Survivin 2B (AYACNTSTL) Monomer is greater than 95% as determined by SEC-HPLC.