

Human HLA-A\*11:01&B2M&KRAS G12C (VVVGACGVGK) Tetramer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM439T-UL


Description	
Source	Recombinant Human HLA-A*11:01&B2M&KRAS G12C (VVVGACGVGK) Tetramer Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus, tetramer is assembled by biotinylated monomer and streptavidin. It contains Gly25-Thr305(HLA-A*11:01), Ile21-Met119(B2M) and VVVGACGVGK peptide.
Accession	AAV53343.1(HLA-A*11:01)&P61769(B2M)&VVVGACGVGK
Molecular Weight	The protein has a predicted MW of 258 kDa. Due to glycosylation, the protein migrates to 260-265 kDa under Non reducing (N) condition based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 90% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Kirsten rat sarcoma 2 viral oncogene homolog (KRAS) is the most commonly mutated oncogene in human cancer. The developments of many cancers depend on sustained expression and signaling of KRAS, which makes KRAS a high-priority therapeutic target. The virtual screening approach to discover novel KRAS inhibitors and synthetic lethality interactors of KRAS are discussed in detail.	

Assay Data

Bis-Tris PAGE



Human HLA-A\*11:01&B2M&KRAS G12C (VVVGACGVGK) Tetramer Protein on Bis-Tris PAGE under Non reducing (N) condition. The purity is greater than 95%.


SEC-HPLC

Human HLA-A\*11:01&B2M&KRAS G12C (VVVGACGVGK) Tetramer Protein, Ultra Low Endotoxin



Cat. No. MHC-HM439T-UL

Assay Data



The purity of Human HLA-A\*11:01&B2M&KRAS G12C (VVVGACGVGK) Tetramer is greater than 90% as determined by SEC-HPLC.