

Biotinylated Human HLA-A*03:01&B2M&gp100 (ALLAVGATK) Monomer Protein



Cat. No. MHC-HM463B

Description

Source	Recombinant Biotinylated Human HLA-A*03:01&B2M&gp100 (ALLAVGATK) Monomer Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Gly25-Thr305 (HLA-A*03:01), Ile21-Met119 (B2M) and ALLAVGATK peptide.
Accession	NP_002107.3(HLA-A*03:01)&P61769(B2M)&ALLAVGATK
Molecular Weight	The protein has a predicted MW of 50.15 kDa. Due to glycosylation, the protein migrates to 51-60 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3-6 months after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Immunization against the ALLAVGATK peptide may occur in melanoma patients and that the immunogenicity of this peptide is similar to that of the majority of melanoma differentiation antigens defined to date. Therefore, the use of ALLAVGATK in active immunotherapy for HLA-A3 melanoma patients, for whom no other immunogenic peptide has been so far described, remains a distinct possibility, providing that patients treated in vivo with this peptide develop a better in vitro response than already demonstrated for HLA-A2 epitopes generated from this same antigenic protein.

Assay Data


Tris-Bis PAGE



Biotinylated Human HLA-A*03:01&B2M&gp100 (ALLAVGATK) Monomer on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Biotinylated Human HLA-A*03:01&B2M&gp100 (ALLAVGATK) Monomer is greater than 95% as determined by SEC-HPLC.