

Biotinylated Human HLA-A*11:01&B2M&HPV16 E7 (IVCPICSQK) Monomer Protein



Cat. No. MHC-HM464B


Description	
Source	Recombinant Biotinylated Human HLA-A*11:01&B2M&HPV16 E7 (IVCPICSQK) Monomer Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Gly25-Thr305 (HLA-A*11:01), Ile21-Met119 (B2M) and IVCPICSQK peptide.
Accession	AAV53343.1(HLA-A*11:01)&P61769(B2M)&IVCPICSQK
Molecular Weight	The protein has a predicted MW of 50.39 kDa. Due to glycosylation, the protein migrates to 52-62 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt.-80°C for 3-6 months after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
HPV16 E7 protein, one of the primary target proteins in molecular targeted therapy for HPV-induced cervical cancer. The affitoxin, ZHPV16E7 affitoxin384 was generated by fusing the modified Pseudomonas Exotoxin A (PE38KDEL) to the HPV16 E7-specific affibody.	

Assay Data


Bis-Tris PAGE



Biotinylated Human HLA-A*11:01&B2M&HPV16 E7 (IVCPICSQK) Monomer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Biotinylated Human HLA-A*11:01&B2M&HPV16 E7 (IVCPICSQK) Monomer is greater than 95% as determined by SEC-HPLC.