

Biotinylated Human HLA-A*02:01&B2M&TdT (ALYDKTKRIFL) Monomer Protein



Cat. No. MHC-HM487B

Description

Source	Recombinant Biotinylated Human HLA-A*02:01&B2M&TdT (ALYDKTKRIFL) Monomer Protein is expressed from HEK293 with His tag and Avi tag at the C-terminus. It contains Gly25-Thr305 (HLA-A*02:01), Ile21-Met119 (B2M) and ALYDKTKRIFL peptide.
Accession	A0A140T913(HLA-A*02:01)&P61769(B2M)&ALYDKTKRIFL
Molecular Weight	The protein has a predicted MW of 50.80 kDa. Due to glycosylation, the protein migrates to 52-65 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	>95% as determined by Bis-Tris PAGE >95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

The terminal deoxynucleotidyl transferase (TdT) belongs to the X family of DNA polymerases. This unusual polymerase catalyzes the template-independent addition of random nucleotides on 3'-overhangs during V(D)J recombination. The biological function and intrinsic biochemical properties of the TdT have spurred the development of numerous oligonucleotide-based tools and methods, especially if combined with modified nucleoside triphosphates.


Assay Data

Bis-Tris PAGE



Biotinylated Human HLA-A*02:01&B2M&TdT (ALYDKTKRIFL) Monomer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Biotinylated Human HLA-A*02:01&B2M&TdT (ALYDKTKRIFL) Monomer is greater than 95% as determined by SEC-HPLC.