

Cat. No. MHC-HM62RT7

Description

Source	Recombinant Alexa Fluor 647-Labeled Human Peptide Ready HLA-B*27:05&B2M Tetramer Protein is expressed from HEK293 with His tag and Avi tag at the C-terminus. It contains Gly25-Pro300 (HLA-B*27:05) and Ile21-Met119 (B2M).
Accession	AAA36221.1(HLA-B*27:05)&P61769(B2M)
Molecular Weight	The protein has a predicted MW of 253.6 kDa.
Wavelength	Excitation Wavelength: 640 nm Emission Wavelength: 672 nm

Formulation and Storage

Formulation	Supplied as 0.22 µm filtered solution in PBS, 100mM L-arginine, 0.2% BSA (pH 7.4).
Storage	Valid for 6 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

HLA-B*27:05&B2M&Peptide ready Monomer is absent from peptide, namely peptide-receptive MHC. It can be loaded with antigenic peptides matching HLA-B*27:05. Peptide ready MHC molecules comprising human HLA alleles and B2M, which can be readily tetramerized and loaded with peptides of choice in a high-throughput manner.