

Mouse c-MPL/Thrombopoietin R Protein

Cat. No. MPL-MM101

Description

Source	Recombinant Mouse c-MPL/Thrombopoietin R Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gln26-Trp482.
Accession	Q08351-1
Molecular Weight	The protein has a predicted MW of 52.3 kDa. Due to glycosylation, the protein migrates to 65-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage

Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

The c-mpl gene encodes a member of the hematopoietic cytokine receptor superfamily. This gene was discovered through the study of a murine retrovirus which induces an acute myeloproliferative syndrome in mice. MPLV (for myeloproliferative leukemia virus) has transduced a truncated and constitutively activated form of the c-mpl receptor chain. High levels of c-mpl expression were detected in 45% of acute myeloid leukemia (AML). No correlation was found between c-mpl expression and the French-American-British classification subtype of AML.


Assay Data

Bis-Tris PAGE



Mouse c-MPL on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Mouse c-MPL is greater than 95% as determined by SEC-HPLC.