

Human MUC18/CD146 Protein

Cat. No. MUC-HM118

Description

Source	Recombinant Human MUC18/CD146 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Val24-Gly559.
Accession	P43121
Molecular Weight	The protein has a predicted MW of 60.8 kDa. Due to glycosylation, the protein migrates to 80-110 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C . Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

METCAM/MUC18 is an integral membrane cell adhesion molecule (CAM) in the Ig-like gene super-family. It can carry out common functions of CAMs which is to perform intercellular interactions and interaction of cell with extracellular matrix in tumor microenvironment, to interact with various signaling pathways and to regulate general behaviors of cells. The above positive correlation is fortuitous; actually METCAM/MUC18 acts as a tumor and metastasis suppressor for the ovarian carcinoma cells.


Assay Data

Bis-Tris PAGE



Human MUC18 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC




The purity of Human MUC18 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Human MUC18, His Tag ELISA

0.1µg Human MUC18, His Tag Per Well



Immobilized Human MUC18, His Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Anti-MUC18, hFc Tag with the EC50 of 12.5ng/ml determined by ELISA (QC Test).