

# Human NGAL/Lipocalin-2 Protein, Ultra Low Endotoxin



Cat. No. NGL-HM201-UL

## Description

<b>Source</b>	Recombinant Human NGAL/Lipocalin-2 Protein is expressed from HEK293 with hFc (IgG1) tag at the C-Terminus. It contains Gln21-Gly198.
<b>Accession</b>	P80188-1
<b>Molecular Weight</b>	The protein has a predicted MW of 47.2 kDa. Due to glycosylation, the protein migrates to 48-52 kDa based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 0.01 EU per µg by the LAL method.
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage


<b>Formulation</b>	Lyophilized from 0.22µm filtered solution in 50mM MES, 150mM NaCl (pH 6.5). Normally 8% trehalose is added as protectant before lyophilization.
<b>Reconstitution</b>	Dissolve the lyophilized protein in 50mM MES (pH 6.5). Please refer to the Certificate of Analysis for detailed instructions.
<b>Storage</b>	-20 to -80°C for 24 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

Acute kidney injury (AKI) is one of the most common complications of various serious conditions, and early diagnosis is therefore critical for the treatment of AKI. Recent evidence demonstrates that neutrophil gelatinase-associated lipocalin (NGAL) is closely associated with AKI. Several experimental and clinical studies have shown that the expression of urine and serum NGAL increases significantly in AKI. NGAL shows potential to be a new effective early biochemical marker of AKI. Further studies are needed to confirm the significant advantages of NGAL in the diagnosis of early AKI and its value in clinical applications.

## Assay Data


### Bis-Tris PAGE



Human NGAL on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC


Assay Data



The purity of Human NGAL is greater than 95% as determined by SEC-HPLC.

ELISA Data

**Human NGAL, hFc Tag ELISA**  
0.1µg Human NGAL, hFc Tag Per Well



Immobilized Human NGAL, hFc Tag at 1 µg/ml (100µl/Well) on the plate. Dose response curve for Biotinylated Anti-NGAL Antibody, hFc Tag with the EC50 21.4ng/ml determined by ELISA.