

Human NKp46/NCR1/CD335 Protein

Cat. No. NKP-HM246

Description

Source	Recombinant Human NKp46/NCR1/CD335 Protein is expressed from HEK293 with hFc tag at the C-Terminus. It contains Gln22-Asn255.
Accession	O76036-1
Molecular Weight	The protein has a predicted MW of 53.2 kDa. Due to glycosylation, the protein migrates to 68-70 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 24 months as supplied from date of receipt. -80°C for 3-6 months after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

NKp46, along with NKp30 and NKp44, are activating receptors that have been collectively termed the natural cytotoxicity receptors (NCR). These receptors lack significant sequence homology to one another. They are expressed almost exclusively by NK cells and play a major role in triggering some of the key lytic activities of NK cells.


Assay Data

Tris-Bis PAGE



Human NKp46 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC




The purity of Human NKp46 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Human NKp46, hFc Tag ELISA

0.05µg Human NKp46, hFc Tag Per Well



Immobilized Human NKp46, hFc Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Anti-Nkp46 Antibody, hFc Tag with the EC50 of 13.1ng/ml determined by ELISA.