

Cynomolgus NPR1/NPRA Protein

Cat. No. NPR-CM301

Description

Source	Recombinant Cynomolgus NPR1/NPRA Protein is expressed from HEK293 with mFc (IgG1) tag at the C-terminus. It contains Gly33-Glu473.
Accession	XP_005541810.1
Molecular Weight	The protein has a predicted MW of 74.56 kDa. Due to glycosylation, the protein migrates to 80-110 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage

Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

NPR1 (natriuretic peptide receptor 1), a receptor of ANP (atrial natriuretic peptide) which acting through NPR1, provokes hypotension. NPR1 was abundantly expressed in endothelial cells and smooth muscle cells of small arteries and arterioles. NPR1 plays a crucial role in ANP-mediated blood pressure regulation, presumably by a mechanism that is RGS2-dependent in the acute phase and RGS2-independent in the chronic phase.


Assay Data

Tris-Bis PAGE



Cynomolgus NPR1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Cynomolgus NPR1 is greater than 95% as determined by SEC-HPLC.