Cynomolgus Neuropilin-1 Protein





Description	
Source	Recombinant Cynomolgus Neuropilin-1 Protein is expressed from HEK293 with His tag at the C-Terminus.
	It contains Phe22-Pro644.
Accession	XP_005564992.1
Molecular Weight	The protein has a predicted MW of 71.01 kDa. Due to glycosylation, the protein migrates to 80-100 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC
F	101

Formulation and Storage


Formulation	Lyophilized from 0.22μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 μ g/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Neuropilin 1 (NRP1) is a transmembrane glycoprotein that acts as a co-receptor for a number of extracellular ligands including class III/IV semaphorins, certain isoforms of vascular endothelial growth factor and transforming growth factor beta.

Assay Data

Tris-Bis PAGE




Cynomolgus Neuropilin-1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Cat. No. NRP-CM101



Assay Data



The purity of Cynomolgus Neuropilin-1 is greater than 95% as determined by SEC-HPLC.