

Cynomolgus PADI4 Protein

Cat. No. PAD-CE104



Description

Source	Recombinant Cynomolgus PADI4 Protein is expressed from E.coli with His tag at the N-Terminus. It contains Met1-Pro663.
Accession	A0A2K5USI5
Molecular Weight	The protein has a predicted MW of 76.6 kDa same as Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Supplied as 0.22µm filtered solution in 20mM Tris, 500mM NaCl, 10% Glycerol, 0.5mM TCEP, 1mM EDTA (pH 8.0).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Peptidylarginine deiminase type4 (PADI4) was firstly identified as a non-MHC RA genetic risk factor. Furthermore, PADI4 risk allele possessed the association with bone damage regardless of anti citrullinated peptide antibody (ACPA) positivity in Asian RA patients. PADI4 gene codes PAD4 protein which has post-translational modification activity (citrullination). Padi4 is mainly expressed in myeloid cells and granulocytes.

Assay Data

Bis-Tris PAGE



Cynomolgus PADI4 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.