

# Biotinylated Mouse uPAR/PLAUR Domain 3 Protein, Ultra Low Endotoxin



Cat. No. PAR-MM1D3B-UL

Description	
Source	Recombinant Biotinylated Mouse uPAR/PLAUR Domain 3 Protein is expressed from HEK293 with His tag and Avi tag at the C-terminus. It contains Phe213-Gly298.
Accession	P35456-1
Molecular Weight	The protein has a predicted MW of 11.83 kDa. Due to glycosylation, the protein migrates to 25-45 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage


Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

The receptor (u-PAR) for urokinase plasminogen activator (u-PA) is a three-domain protein, GPI-anchored to the cell surface, which focuses the enzymatic activity of u-PA, and allows the cell surface activation of plasminogen. Regulation of the activity of u-PA is also mediated by u-PAR.

## Assay Data


### Bis-Tris PAGE



Biotinylated Mouse uPAR Domain 3 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC

Assay Data



The purity of Biotinylated Mouse uPAR Domain 3 is greater than 95% as determined by SEC-HPLC.