

Mouse PSGL-1 Protein

Cat. No. PSG-MM162

Description

Source	Recombinant Mouse PSGL-1 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gln42-Cys307.
Accession	Q62170
Molecular Weight	The protein has a predicted MW of 28.74 kDa. Due to glycosylation, the protein migrates to 68-108 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

P-selectin glycoprotein ligand-1 (PSGL-1) has long been studied as an adhesion molecule involved in immune cell trafficking and is recognized as a regulator of many facets of immune responses by myeloid cells. PSGL-1 also regulates T cell migration during homeostasis and inflammatory settings.

Assay Data

Bis-Tris PAGE



Mouse PSGL-1 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.