

Human PSMA/FOLH1 Protein

Cat. No. PSM-HM110

Description

Source	Recombinant Human PSMA/FOLH1 Protein is expressed from Expi293 with His tag at the N-terminal. It contains Lys44-Ala750.
Accession	Q04609-1
Molecular Weight	The protein has a predicted MW of 80.6 kDa. Due to glycosylation, the protein migrates to 90-110 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 $\mu\text{g}/\text{ml}$ is recommended (usually we use 1mg/ml solution for lyophilization). Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -20 to -80°C for 3-6 months in unopened state after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please avoid freeze-thaw cycles.

Background

Prostate-specific membrane antigen (PSMA) is an enzyme that in humans is encoded by the FOLH1 (folate hydrolase 1) gene, also known as Glutamate carboxypeptidase II (GCPII). Human PSMA is highly expressed in the prostate, roughly a hundred times greater than in most other tissues. In some prostate cancers, PSMA is the second-most upregulated gene product, with an 8- to 12-fold increase over levels in noncancerous prostate cells.

Assay Data


Tris-Bis PAGE



Human PSMA on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Human PSMA is greater than 95% as determined by SEC-HPLC.

ELISA Data

Human PSMA, His Tag ELISA

0.2µg Human PSMA, His Tag Per Well



Immobilized Human PSMA, His Tag at 2µg/ml (100µl/Well) on the plate. Dose response curve for Anti-PSMA Antibody, hFc Tag with the EC50 of 15.2ng/ml determined by ELISA.