

Human SEZ6 Sushi4 Domain Protein



Cat. No. SEZ-HM20D

Description	
Source	Recombinant Human SEZ6 Sushi4 Domain Protein is expressed from HEK293 with hFc tag at the C-terminus. It contains Gln769-Leu832.
Accession	Q53EL9-1
Molecular Weight	The protein has a predicted MW of 32.94 kDa. Due to glycosylation, the protein migrates to 40-48 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Seizure-related protein 6 (Sez6) contributes to chronic pain development as sez6 knockout mice show attenuated pain behaviours after peripheral nerve injury, compared with control mice. The type I transmembrane isoform of Sez6 is cleaved by the β-amyloid precursor protein cleavage enzyme 1 (BACE1), resulting in Sez6 extracellular domain shedding from the neuron surface.	

Assay Data

Bis-Tris PAGE



Human SEZ6 Sushi4 Domain on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Human SEZ6 Sushi4 Domain Protein

Cat. No. SEZ-HM20D



Assay Data



The purity of Human SEZ6 Sushi4 Domain is greater than 95% as determined by SEC-HPLC.