

Human SFRP1/SARP2 Protein

Cat. No. SFP-HM101

Description

Source	Recombinant Human SFRP1/SARP2 Protein is expressed from Expi293 with His tag at the C-terminal. It contains Ser32-Lys314.
Accession	Q8N474
Molecular Weight	The protein has a predicted MW of 33.63 kDa. Due to glycosylation, the protein migrates to 43-48 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE

Formulation and Storage

Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4). Please dilute to the desired concentration according to the concentration of the solution shown on the product label.
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please do not repeated freeze-thaw cycles.

Background

Secreted frizzled-related protein 1 (SFRP1) is a gene that belongs to the secreted glycoprotein SFRP family. SFRP1 has been classified as a tumor suppressor gene due to the loss of expression in various human cancers, which is mainly attributed by epigenetic inactivation via DNA methylation or transcriptional silencing by microRNAs.

Assay Data

Tris-Bis PAGE



Human SFRP1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.