

Biotinylated Human SERPINF2/A2AP Protein

Cat. No. SPF-HM401B

Description

Source	Recombinant Biotinylated Human SERPINF2/A2AP Protein is expressed from HEK293 with His tag and Avi tag at the C-terminus. It contains Met28-Lys491.
Accession	P08697-1
Molecular Weight	The protein has a predicted MW of 54.62 kDa. Due to glycosylation, the protein migrates to 60-75 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Large-pool solvent/detergent (SD) plasma for transfusion exhibits reduced alpha 2-antiplasmin (alpha2-AP; SERPINF2) functional activity. The reason for the loss of alpha2-AP has not been described and could be due to the SD incubation itself and/or to the processing steps implemented to remove the solvent and the detergent.


Assay Data

Bis-Tris PAGE



Biotinylated Human SERPINF2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC




The purity of Biotinylated Human SERPINF2 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Biotinylated Human SERPINF2, His Avi Tag ELISA

0.05µg Biotinylated Human SERPINF2, His Avi Tag Per Well



Immobilized Biotinylated Human SERPINF2, His Avi Tag at 0.5µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Anti-SERPINF2 Antibody, hFc Tag with the EC50 of 13.1ng/ml determined by ELISA.