

Biotinylated Human Mature TGF beta 2 Protein



Cat. No. TG2-HM00MB

Description	
Source	Recombinant Biotinylated Human Mature TGF beta 2 Protein is expressed from HEK293 with Avi tag at the C-Terminus. It contains Ala303-Ser414.
Accession	P61812-1
Molecular Weight	The protein has a predicted MW of 14.53 kDa same as Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage	
Formulation	Supplied as 0.22µm filtered solution in 4mM HCl.
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Transforming growth factor beta(2) (TGF-beta(2)), a growth regulator of human lens epithelial cells (HLECs), also regulates the death of these cells. TGF-beta(2)-induced apoptosis in HLECs was preceded by an induction of reactive oxygen species (ROS) and a decrease in glutathione in the intracellular content, indicating that this factor induces oxidative stress in HLECs.	

Assay Data

Bis-Tris PAGE




Biotinylated Human Mature TGF beta 2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

Assay Data

Biotinylated Human Mature TGF beta 2, Avi Tag ELISA

0.5µg Biotinylated Human Mature TGF beta 2, Avi Tag Per Well



Immobilized Biotinylated Human Mature TGF beta 2, Avi Tag at 5µg/ml (100µl/well) on the streptavidin precoated plate (5µg/ml). Dose response curve for Human TGF-beta RII, mFc Tag with the EC50 of 0.17µg/ml determined by ELISA (QC Test).