

Mouse TIGIT Protein, Ultra Low Endotoxin

Cat. No. TIG-MM210-UL

Description

Source	Recombinant Mouse TIGIT Protein is expressed from HEK293 with hFc (IgG1) tag at the C-Terminus. It contains Gly26-Thr143.
Accession	NP_001139797.1
Molecular Weight	The protein has a predicted MW of 39.7 kDa. Due to glycosylation, the protein migrates to 48-57 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

TIGIT, also called Vstm3, Vsig9, and WUCAM, is a transmembrane protein in the CD28 family of the Ig superfamily proteins. IGIT is expressed at low levels on peripheral memory and regulatory CD4 T-cells and NK cells and is up-regulated following activation of these cells.

Assay Data


Bis-Tris PAGE



Mouse TIGIT on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Mouse TIGIT is greater than 95% as determined by SEC-HPLC.

ELISA Data

Mouse TIGIT, hFc Tag ELISA

0.5µg Mouse CD155, His Tag Per Well



Immobilized Mouse CD155, His Tag at 5µg/ml (100µl/Well) on the plate. Dose response curve for Mouse TIGIT, hFc Tag with the EC50 of 56.5ng/ml determined by ELISA