

Human TIM-3/HAVCR2 Protein, Ultra Low Endotoxin



Cat. No. TIM-HM131-UL


Description	
Source	Recombinant Human TIM-3/HAVCR2 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Ser22-Arg200.
Accession	Q8TDQ0-1
Molecular Weight	The protein has a predicted MW of 20.7 kDa, Due to glycosylation, the protein migrates to 50-60 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.001 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt.-80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Hepatitis A virus cellular receptor 2 (HAVCR2), also known as T-cell immunoglobulin and mucin-domain containing-3 (TIM-3), is a protein that in humans is encoded by the HAVCR2 gene.TIM3 is an immune checkpoint and together with other inhibitory receptors including programmed cell death protein 1 (PD-1) and lymphocyte activation gene 3 protein (LAG3) mediate the CD8 T-cell exhaustion. TIM3 has also been shown as a CD4 Th1-specific cell surface protein that regulates macrophage activation and enhances the severity of experimental autoimmune encephalomyelitis in mice.	

Assay Data


Bis-Tris PAGE



Human TIM-3 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Human TIM-3 is greater than 95% as determined by SEC-HPLC.

ELISA Data

Human TIM-3, His Tag ELISA

0.05µg Human TIM-3, His Tag Per Well



Immobilized Human TIM-3, His Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Anti-TIM-3 Antibody, hFc Tag with the EC50 of 4.3ng/ml determined by ELISA.