

Mouse TIMP-1 Protein

Cat. No. TIM-MM201

Description

Source	Recombinant Mouse TIMP-1 Protein is expressed from HEK293 with hFc tag at the N-Terminus. It contains Cys25-Arg205.
Accession	P12032
Molecular Weight	The protein has a predicted MW of 46.38 kDa. Due to glycosylation, the protein migrates to 55-65 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Tissue inhibitor of metalloprotease-1 (TIMP-1) is a tissue inhibitor of matrix metalloproteinases (MMPs). It however exerts multiple effects on biological processes, such as cell growth, proliferation, differentiation and apoptosis, in an MMP-independent manner.


Assay Data

Bis-Tris PAGE



Mouse TIMP-1 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Mouse TIMP-1 is greater than 95% as determined by SEC-HPLC.