

Mouse VLDLR Protein

Cat. No. VDR-MM101



Description	
Source	Recombinant Mouse VLDLR Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gly28-Ser797.
Accession	P98156
Molecular Weight	The protein has a predicted MW of 86.1 kDa. Due to glycosylation, the protein migrates to 110-130 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage	
Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
VLDLR cerebellar hypoplasia (VLDLR-CH) is characterized by non-progressive congenital ataxia that is predominantly truncal and results in delayed ambulation, moderate-to-profound intellectual disability, dysarthria, strabismus, and seizures.VLDLR-CH is inherited in an autosomal recessive manner. Carrier testing for at-risk relatives, prenatal testing for a pregnancy at increased risk and preimplantation genetic testing are possible when the pathogenic variants in a family are known.	


Assay Data

Bis-Tris PAGE



Mouse VLDLR on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SPR Data



Mouse VLDLR, His Tag immobilized on CM5 Chip can bind Mouse PCSK9, His Tag with an affinity constant of 0.28 nM as determined in SPR assay (Biacore T200).