

# Mouse VEGF R2/KDR Protein

Cat. No. VGF-MM1R2

## Description

<b>Source</b>	Recombinant Mouse VEGF R2/KDR Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Ala20-Glu762.
<b>Accession</b>	P35918-1
<b>Molecular Weight</b>	The protein has a predicted MW of 84.11 kDa. Due to glycosylation, the protein migrates to 110-140 kDa based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 1EU per $\mu\text{g}$ by the LAL method.
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage

<b>Formulation</b>	Lyophilized from 0.22 $\mu\text{m}$ filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
<b>Reconstitution</b>	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 $\mu\text{g}/\text{ml}$ is recommended. Dissolve the lyophilized protein in distilled water.
<b>Storage</b>	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

Tyrosine-protein kinase that acts as a cell-surface receptor for VEGFA, VEGFB and PGF, and plays an essential role in the development of embryonic vasculature, the regulation of angiogenesis, cell survival, cell migration, macrophage function, chemotaxis, and cancer cell invasion. The tyrosine kinase receptor vascular endothelial growth factor receptor 2 (VEGFR2) is a key regulator of angiogenesis.


## Assay Data

### Bis-Tris PAGE



Mouse VEGF R2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC




The purity of Mouse VEGF R2 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Mouse VEGFR2, His Tag ELISA

0.5µg Human VEGF165, No Tag Per Well



Immobilized Human VEGF165, No Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Mouse VEGFR2, His Tag with the EC50 of 1.15µg/ml determined by ELISA.